

DNA Isolation Tools of Hain Lifescience

GenoType[®] DNA Isolation Kit

Hain Lifescience's **GenoType[®] DNA Isolation Kit** allows the quick isolation of highly purified DNA. The yield DNA is suitable for any usual downstream amplification reaction (PCR) from smallest amounts of human whole blood, body fluids, tissue samples or eukaryotic cell cultures.



Your benefits at a glance:

- DNA isolation and purification in approximately 12 minutes
- Reliable DNA extraction from smallest amounts of source material (1-50 μ l)
- Highest yields (1-2 μ g DNA)
- Convenient handling – only a few steps for quick and optimum results

GENO•CARD®

Hain Lifescience's **GENO•CARD®** allows the quickest DNA isolation ever. One drop of the sample to be investigated is spotted on the **GENO•CARD®**. Subsequently, a small piece of the carrier matrix is punched out and directly used in the amplification reaction without further purification. **GENO•CARD®** provides a simple and time-saving way of DNA isolation.



Your benefits at a glance:

- DNA isolation of amplifiable DNA in only a few minutes
- Simplest handling
- Storage of spotted material for months

Hain Lifescience GmbH

Hardwiesenstraße 1 | 72147 Nehren | Germany
Tel.: +49 (0) 74 73- 94 51- 0 | Fax: +49 (0) 74 73- 94 51- 99
E-Mail: info@hain-lifescience.de | www.hain-lifescience.de

